

SLIDE HR®



V-Shape
PLATINUM



COOL GLAMOR®



PZLASER®



EosICE®
MAX



10D
COOL GLAMOR®



MDSculpt iD™

EosICE®
Pro Max



EosICE®
PLATINA



EosICE™
P
R
O



ThemaLift RF®



UMS Sculptor®



EosICE®

PicoIris®



Fractional Master™



UMSELLA™



Microkits.cl
Servicios Tecnológicos

Mauricio Moreno

Email: mmoreno@microkits.cl

Cell/WhatsApp: +56 9 62681766

Innovation
Never Stop
All to create value
for **you**

Let the World Love

MEDICAL



TGA

ISO 13485

UMSFACE

To Make the Best Skin of Your Life





Treatment Areas

- ◆ Forehead
- ◆ Eyes
- ◆ Eyebrows Cheeks
- ◆ Jawline
- ◆ Chin
- ◆ Cervical Spine

Technology Behind UMS FACE


UMS FACE is the first device on the market utilizing the simultaneous application of both the synchronized RF and HIFES energies. Synchronized RF remodels & smooths skin by heating the dermis and increasing the levels of collagen and elastin fibers.

While the HIFES technology restores & elevates support of facial tissues by selectively contracting muscles and increasing density & quality of muscle structure.



KEY BENEFITS For UMS FACE

- ◆ **2 THERAPIES IN ONE PROCEDURE**
Combination of synchronized RF and HIFES energies for simultaneous treatment of facial skin and muscles.
- ◆ **NEEDLE-FREE**
The procedure is completely non-invasive without any downtime.
- ◆ **HANDS FREE**
The procedure runs independently once applicators are affixed.
- ◆ **FULL FACE IN 20 MINUTES**
UMS FACE procedure is comfortable and operator-friendly treatment lasting for only 20 minutes.



V-Shape PLATINUM


Make Yourself Smoother and Sexier From Inside to Outside

A New Star Devised and Crafted by PZ LASER


Application

- Fat reduction
- Health regimen
- Texture improvement
- Wrinkle removal
- Cellulite reduction

• 5 ACTIONS IN ONLY ONE TREATMENT



Remodeling action
Reshapes the silhouette and smooths out wrinkles.



Muscle action
Tones muscles and firms' tissues.



Magic Ball Applicator

It uses an innovative COMPRESSIVE MICROVIBRATION technology to feel the body and produce a customized treatment to tone, tighten and smooth out imperfections.



Large Applicator

Four in 1 technology works with a combination of infrared light, bipolar radio frequency energy, mechanical rollers and vacuum.

- Vacuum + Mechanical Rollers
- Bipolar RF
- Infrared Light



Medium Applicator

Four in 1 technology and 1 hand-pieces with 2 changeable applicators.

The medium handle have two different heads , combine with RF , Negative pressure, Vacuum infrared light , rollers technology together.



Creative Small Applicator

One hand-piece with 3 different sizes for different treatment areas, easy to operate.

A well-deserved facial carving artifact.



EosICE[®]

Pro Max

NEWEST REVOLUTIONARY HAIR REMOVAL DIODE LASER DEBUTS GLOBALLY

- Germany LIMO Collimator Lenses
- 360° Cooling Peltier
- Germany LIMO Focusing Lens




Using FocusLight assembled Germany Dilas Laser bar

- The effect is increased by 200%
- The pain is reduced by 70%

Android Intelligent UI

 15.6 inch Luxury Android Smart Screen
Remote control, unlimited upgrades and storage



New
DEBUTS
GLOBAL
LY

AI Diode Laser Technology

The world's first AI diode laser technology, equipped with hair detection to intelligently recommend the best parameters.




AI Hair Analyzer



3 Changeable Spot Sizes Handle

① One handle three spot

15*15mm, 15*30mm 8*8mm can be switched freely



② Four in1 Personalized Platform

808nm, 755nm, 940nm, 1064nm wavelengths, from single wave to 4 bands, what you want is what I will do.



World 's First Optical Waveguid Nose Tip

First Nose tips with optical waveguide jointly developed by German LIMO

+100% Treatment Gentleness **+100%** Treatment Effect



EosICE[®] MAX

SMARTER, STRONGER, FASTER

- ✓ Alex 755nm
- ✓ Diode 808nm
- ✓ Infrared 940nm
- ✓ YAG 1064nm

Technology Advantages



Exclusive Patent

Four-Wavelength Laser, multi-layer energy penetration achieve the purpose of hair removal, clean and comfortable.



Germany Laser Bar

The laser transmitter imported from the USA, and the service life is not less than 20 million times.



First Annual Cooling Pettier Insides Handle

The American original high-power semiconductor electronic cooler, water pump imported from Italian provides continuous cooling for the laser to ensure operation for 24 hours.



-18°C Cooling

Adopting TEC cooling system, contact cooling, comfort enhancement.



15.6 Inch Smart Android Operation system

The intelligent system automatically recommends appropriate parameters according to skin color, hair softness and hair color, which is quick and effective.





Diameter 8 mm



15*30mm



15*15mm



1.5 times work efficiency: 72 hours of super long work to ensure stable function in busy clinic/salons.



Eos*ICE*[®] PLATINA

4 in 1 PERSONALIZED PLATFORM




Innovative 4D Tech Laser Hair Removal

808nm, 755nm, 940nm, 1064nm wavelength

From single wavelength to 4 wavelengths, we can customize whatever you want.



Precisely target melanin in hair follicles, inhibit hair regeneration in different parts of the human body



In orderto provide you with optimal treatment possibilities, the portable Eos ICE Platina comes with:

- The "Gold Standard" 4D Tech for Hair Removal
- Unique Built-in 4 in 1 cooling
- Up to 15Hz Repetition - In Motion Working Mode
- Intuitive software with preset parameters and professional mode

The unique features of the Eos ICE Platina, enable you to provide your patients with:

- Short treatment time
- Maximal treatment comfort
- Optimal results
- Skin type (I-VI) & Hair color and texture adjustable

Ultra PicoIris[®]

Highest Power | Shortest Pulse Width | Versatility for Difficult Cases



Versatility for Difficult Cases



Tattoos
Removal



Freckle
Removal



Pigment
Removal



Hair
Removal



Veins



Onychomycosis



Technology Behind PicoIris

PicoIris is the next-generation product following the top-of-the-line series of Q-switched Nd:YAG laser devices perfected by PZ LASER, challenging the pulse mode from nanosecond to picosecond.

PicoIris offers 5 wavelengths: 532nm, 1640nm, 755nm, 585nm, 650nm and Fractional handle, which gives the versatility to treat even the most difficult cases.



Break-through 5 working modes

-  **PICOSECOND**
Pigments, Tattoos, Birthmarks, Skin Rejuvenation
-  **PTP**
Melasma, Moles, Tattoo, Hyperpigmentation.
-  **MULTI PULSE**
Chloasma, Stubborn Pigments.
-  **LONG PULSE**
Vascular, Scars, Acne, Rosacea, Diffuse Redness, Onychomycosis, Rejuvenation, Firming, Whitening, Shrinking Pores, Remove Fine Wrinkles, Hair Removal, Improve Veins.
-  **LONG PULSE PRO**
Vascular, Scars, Acne, Rosacea, Diffuse Redness, Onychomycosis, Rejuvenation, Firming, Whitening, Shrinking Pores, Remove Fine Wrinkles, Hair Removal, Improve Veins.



360° COOLING A SCULPT YOU WITH GLAMOR!



Features

-  Patented Massage Mode
-  Patented Pre-heating Lymph Mode
-  Constant Temperature Technology
-  For ALL Body Profiles and Concerns
-  Diverse Range of Applicators



REINVENTS FAT FREEZING MACHINE

COOL GLAMOR uses state-of-the-art freezing technology removes excess subcutaneous fat from a specific location. Designed with 360° Surround Cooling Technology, the multiple cooling cups can fit all curves and contours from as high as your chin all the way down to the knees, to deliver an absolute non-invasive fat reduction way.

Treat 4 Areas Simultaneously in Only 40 Minutes

- ◆ COOL GLAMOR will come with 4 handles can work at the same time
- ◆ Each handle is individually controlled and can be adjusted with different parameters.
- ◆ Available to treat multiple areas or multiple people simultaneously

Thoughtful Fast and Easy Applicator Exchange



Make a 360 Sculpt You with Glamor



COOL GLAMOR[®] II

THE WORLD'S LEADING CRYOLIPILYSIS FAT FREEZING MACHINE

8 innovative designs

- 7 in 1 multi-size applicator
- Medical flexible LSR silicone cups
- innovative "Preheat" technology
- Patented massage mode
- Innovative constant temperature technology
- Upgrade vacuum system, the suction is increased by 20%
- Self-cleaning
- 2 applicatos could work at the same time or individually



Treat 2 Areas Simultaneously in Only 40 Minutes

The COOL GLAMOR is equipped with an expandable range of applicators that makes it easier than ever to achieve significant results on the most challenging of areas.

Upgraded Technology for Enhanced Results.

The upgraded 360° surround cooling of the COOL GLAMOR delivers cooling energy throughout the entire applicator and makes it at a faster rate, thus increasing efficiency by up to 25%.

Humanized treatment advantage

-  Patented Massage Mode
-  Patented Pre-heating Lymph Mode
-  Constant Temperature Technology
-  For ALL body profiles and concerns
-  Diverse Range of Applicators




COOL GLAMOR[®] pro

ELECTROMAGNETIC+ FREEZING

HIFEM + 360° CRYOLIPOLYSIS

-  **Patented technology,** high strength, high focus.
-  **Double effect,** fat-dissolving & muscle-building in one step.
-  **Non-invasive,** safe and comfortable without recovery period.
-  **Multiple parts work at one time,** saving time and convenience.
-  **After treatment,** the effect of muscle increase continued to 19%.




Cooling 360° Fat Freezing



Highly Focused Electromagnetic Build Muscle



Combining two technologies



Why choose us ?

The first in the industry: dual-energy working system (refrigeration technology + high focus electromagnetic technology)

Powerful energy: dual power supply configuration, independent energy output is stable, no attenuation.



UMS Sculptor®

Non-invasive Fat Dissolving
Muscle Building
Private Anti-Age

Technology Upgrade

150HZ, more than 20,000 "ultra-muscle exercises", achieving extraordinary performance.

Function Upgrade

For muscle increase, rectus abdominis and pelvic floor muscle recovery.

Configuration Upgrade

The 4 applicators have independent intelligent working modes. No water cooling requested, offering convenient maintenance.

Design Upgrade

The shell adopts high-level aesthetic color and the shape is outstanding.

Ergonomic Upgrade

Intelligent full display and the applicators have a mechanics design.



Multi-Function, New Technology

① (Fat Burning and Muscle Building) Innovative Hifem Technology

② (Private Health) Patented Electromagnetic Technology




SELECTIVE PURCHASE

SEATED WORKING HEAD

High-energy-focused electromagnetic technology fully activates pelvic floor muscle. It perfectly simulates the curve of the buttocks and accurately targets the gluteus. The 15° inclination design allows the legs to relax naturally.



Nicely Fit in Multiple Mode's



Advantage

1. Fat Reduction+Muscle Building + Private Recovery, satisfying all need of consumption.
2. Non-Surgical, Non-Invasive, No-Downtime
3. Intelligent interface, convenient for operating and using.
4. Original German imported, the effect and quality are guaranteed.
5. The clinically proven effect is obvious and lasting.


UMSELLA™

Non invasive | Private Anti-Aging

Principle

High-intensity focused electromagnetic waves stimulate the pelvic floor muscle

Electromagnetic Energy Focused→Pulse Deep Penetration→Stimulate Pelvic Floor Muscles→1 Treatment→12,000 Pelvic Floor Muscle Contractions



- **Pelvic Muscle Rejuvenation**
- **Comfortable Treatment**
- **Obvious Effect**

Dual Mechanism


- **Sacral Nerve Stimulation:** Through the high-intensity magnetic field, it directly acts on the sacral nerve, promotes the repeated synchronization of neurons before and after the synapse, strengthens the nerve muscle reflex arc, and improves the nerve's innervation of the target organ.
- **Pelvic Muscle Stimulation:** The pulsed electromagnetic field penetrates muscles and fat tissues without attenuation, effectively stimulates the peripheral motor nerves of the pelvic floor, triggers nerve depolarization, generates induced currents and transmits them to the motor endplates, triggers repeated contractions of muscles, thereby strengthening the pelvic floor muscle strength and muscle endurance.

Application of UMSella

- **Comprehensive improvement**
Pelvic floor dysfunction (PFD), stress urinary incontinence, pubococcygeous enhancement
- **Post-partum Recovery**
Improve postpartum problems, metroptosis, quality of life improvement, urine leakage etc.
- **Womem's Health**
Rehabilitation of pelvis muscles, women's private health improvement

Advantage of UMSella

- **Ergonomic Design:** The arc-shape fits in the waist-hip curve and achieves the effect directly.
- **Electromagnetic Stimulation:** Stimulate deeper and wider pelvic floor muscles and nerves.
- **24/7 Operation:** 72- hours work without stopping, high efficiency.
- **Customized Data System:** Make your own personalized treatment plan based on clinical needs.



Quality of life improvement



Prevent stress urinary incontinence



Post-partum Recovery



Rehabilitation of Pelvic Muscles



Pubococcygeous Enhancement



Women's Private health Improvement



A Platform for Body and Face



MRF-X

MRF-X is an advanced workstation for both face and body treatment based on STCT (Smart Temperature Controlling Technology)



Technology behind MRF-X

MRF-X uses RF energy waves with STCT (Smart Temperature Controlling Technology) to heat the deep layer of your skin known as your dermis. This heat stimulates the production of collagen, collagen is the most common protein in your body, so MRF-X will create the framework of your skin and gives your skin its firmness.

Face
Applicator



Body
Applicator



Advantages

- ◆ Smart Temperature Controlling Technology
- ◆ High power 500W
- ◆ User Friendly UI
- ◆ Easier Operation
- ◆ Immediate Results

Appication

- ◆ Cellulite Reduction
- ◆ Body Slimming
- ◆ Weight Loss
- ◆ Skin Tightening
- ◆ Body Shaping
- ◆ Wrinkles Removal
- ◆ Face Lift

HydraFacial™

Get the Best Skin
of Your Life



Our non-invasive treatment improves skin health, addresses individual concerns, and creates a glow like no other. Unlike traditional facial skin care treatments, it's performed via device rather than by hand, for deeper penetration and longer lasting results.

Versatile Applicators - 10 in 1 Platform



•AQUA PEEL



•OXYGEN JET



•EMS FACE



•COLD-HEAT



•MICRO LIFT



•OXYGEN




•ULTRASOUND



•RF




•CLEANSING
SHOVEL



•ANALYZER

10D COOL GLAMOR[®]


To make the magic happen with
your face and body



Various cartridges up to 10 types

Customized treatment depending on the skin type, treatment area and depth of dermis to subcutaneous fat layer.

- ◆ FL(Focused Ultrasound Linear)6 types : Face 4 types (1.5mm+2.0mm+3.0mm+4.5mm with 7MHZ), Body 2 types (6.0mm+9.0mm with 4MHZ).
- ◆ SD(Synergy Dotting)4 types : Face 4 types(1.5mm+2.0mm +3.0mm+4.5mm, can adjust the penetration depth manually).



Patented safety of Coolglamor

FL contact sensor: The cartridge will stop working automatically to avoid damage if not contact with the skin correctly.

FL movement sensor: The movement sensor prevents consecutive shots in one place.

SD gyro sensor: Use SD gyro sensor to prevent overlapping on the same spot and move forward.



Benefits from Coolglamor

- ◆ Provide new experience and effect.
The synergy of MFU provide better results.
- ◆ Profit making device with new MFU technology and various handpieces.
High treatment price can be formed based on the synergy of MFU.
Not only skin but also body contouring is possible with the specialized body applicators.
- ◆ Shorten procedure time.
Improved treatment speed with various of applicators is no dead angle in operation.

Fractional *Master*TM

The Most Powerful Co2 Fractional Laser
for Aesthetic Treatment



Co2 Fractional Laser

- 

Skin Rejuvenation
& Tightening
- 

Vaginal
Tighten
- 

Scar
Recovery
- 

Stretch
Marks
- 

Pigmentation



Co2 Fractional Laser

The fractional CO2 laser takes the CO2 laser beam and fractionates, or pixellates that beam into thousands of tiny little shafts of light. These tiny shafts of light penetrate into the deeper layers of the skin. The skin then repairs those tiny shafts by pushing out the old sun damaged skin and replacing it with new skin.

Treatment Features

① Non-Ablative SkinTight Applicator



② Vulva Mode



④ Fractional Mode

③ Super Pulse
Mode



Highly intelligent analysis of freckle
and acne marks, stretch marks, skin
rejuvenation, and scar repair.



High Configurations

- 

Lasers imported from the USA with sufficient and stable energy
- 

Galvo mirror imported from Germany with balanced micro-beam energy
- 

7-joint arm imported from South Korea, flexible operation and precise output

Perfect LightTM

SAFE, FAST, EFFECTIVE

Platform with two innovative technologies: DPL Technology + NIR Technology Multi-functional choice for both hair removal and skin treatments



APPLICATION

Multi-functional choice for both hair removal and skin treatments, such as pigment removal, freckles removal, vascular removal, acnes removal and skin rejuvenation, etc.


HDPL for Hair

Wavelength 640-850nm

Fast hair removal adopts 640nm-850nm wavelength, a very precise narrow band spectrum for accurate and pain-free hair treatment. It is a technology as effective as laser treatment but with the safety, versatility and cost of pulsed light.

NIR for Face and Body

Wavelength 900-1800nm



NIR uses a high power light source in the near infrared spectrum to achieve safe and effective sub-dermal heating, stimulating the growth of new collagen and tightening lax skin while also imparting a youthful glow to the skin. It can be used on the face or focused on parts of the body. Procedures with the NIR module may be performed via a stationary delivery method or using the in-motion technique for gradual heat buildup and comfortable treatment.



• Body



• Face

SDPL for Skin

Wavelength 530-850nm

Targets all layers of hyper-pigmentation Including dark skin, for all types of pigmented and vascular lesions.

Additional Choice

QBeam^{Plus}TM

The Most Comprehensive Q-Switch Nd:YAG
Laser Treatment Solutions



SUPERIORITY

The Most comprehensive Q-Switch Nd:YAG Laser Treatment solutions available today.

A Single laser module treats 3 separate indications for all skin types, with new and unique technological innovations for each indication.

Technology Behind Q-Beam Plus

Patented “One Lamp, Two Rods” Q-switched Nd: Yag laser mode, which radiates the light of specific wavelength in high peak energy pulses. In 6ns the light penetrates the tissue and is absorbed by the pigmentation resulting in an instantaneous blast.

High Configuration of Q-Beam

- ◆ 3 Kinds of treatment tips are available for our patented Q+ treatment handle, provides you more skillful and easy operations.
- ◆ High power leads to high energy density. Which enables perfect result of endogenous tattoos and pigments removal.



532nm Tip



1064nm Tip



1320nm Tip

Versatility of Treatments

◆ Tattoo Removal

Black, Blue, Brown, red different kinds of old and new tattoos

◆ Pigmented Lesions in Epidermis and Dermis

Age spot, Coffee spot, Sunspot, Freckle, Melasma, Ota, Zygomaticus, Nevus, Mongolian Spot, Pigmentation after inflammation

◆ PTP Mode

Skin rejuvenation, Skin whitening, Skin tightening, Skin texture toning.

Platinum Skin™

A Hair Removal Device that can
Rejuvenate the Skin



Superiority

- ◆ Thoughtful and Innovative Tech
- ◆ A Hair Removal Device that can Rejuvenate the Skin
- ◆ Milk Light + Diode Laser
- ◆ Powerful Choice for both Hair Removal and Skin Rejuvenation

Technology of Platinum Skin


Thoughtful and Innovative Tech with Milk Light + Diode Laser, it is Powerful Choice for both Hair Removal and Skin Rejuvenation

4D Diode Laser Technology

4D Tech Diode Laser offers breakthrough hair removal technology using a diode laser based on an innovative 4 wavelengths. The powerful absorption that is possible with the 808nm+755nm+940nm+1064nm wavelength combined with the contact cooling, treatment coverage, comfort and low maintenance of the diode laser provides an ideal solution for the widest range of hair types and colors.




Precisely target melanin in hair follicles, inhibit hair regeneration in different parts of the human body



Milk Light Technology

Milk Light (NIR) uses a high power light source in the near infrared spectrum to achieve safe and effective sub-dermal heating, stimulating the growth of new collagen and tightening lax skin while also imparting a youthful glow to the skin. It can be used on the face or focused on parts of the body.



Eos Ice Prime

A breakthrough hair removal platform



Applicator
features




Lighter
operation



Touch
cooling system



Intelligent
configuration



• One handle double spot to switch freely.

• Magnetic spot replacement port, seamless fit and easy operation.

• The latest generation of laser technology, ultra Light, small and strong energy.

For all skin types

Effective results for light or dark skin as well as coarse or fine hair.



Intelligent parameter system

The professional data system automatically adjusts parameters according to different needs, different skin colors, genders, parts and the color, thickness and area of the hair to achieve perfect control.

Application

Complete and effective solution for all skin tones, hair types and target areas (Body & Face)

Face: chin and lips, beard.

Body: armpit, arms, leg hair, chest hair, bikini area.



Ultra PicoIris[®]

Photocoustic physical blast wave

A real sense of ultra light picosecond

Application

PICOSECOND

It can be applied to most of the spots, birthmarks or tattoos, skin rejuvenation treatment, tattoo!

PTP

Melasma, Moles, Tattoo, Dark Eye Circles, Pores, Hyperpigmentation.

MULTI PULSE

Mainly used for chloasma, Stubborn stain and pigments.

LONG PULSE

Vascular, Scars, Acne, Rosacea/diffuse redness, Onychomycosis, Rejuvenation, firming, whitening, shrinking pores, Remove fine wrinkles, Heating dermal cells, hair removal, improve Varicose veins.

LONG PULSE PRO

Vascular, Scars, Acne, Rosacea/diffuse redness, Onychomycosis, Rejuvenation, firming, whitening, shrinking pores, Remove fine wrinkles, Heating dermal cells, hair removal, improve Varicose veins.



Hair removal



Improve varicose veins



Onychomycosis



Freckle removal



Acne clearance



Facial rejuvenation



Tattoos removal

Superiority of PicoIris

PicoIris is the next-generation product following the top-of-the-line series of Q-switched Nd:YAG laser devices perfected by PZ LASER over the past 18 years, and challenging the pulse mode from nanoseconds to picoseconds. PicoIris offered 5 wavelengths: 532 nm, 1640nm, 755nm, 585 nm and 650 nm, This variety of options gave it the versatility to treat even the most difficult cases.



Uniqueness of PicoIris

01 Five treatment modes

Innovative 5 treatment modes, make versatility for more difficult cases.

02 User Friendly UI with clinical database

Clinical database recommendations as reference, making it the easiest for even new beginners to start.

03 Integrated power supply and 4 operate independently control systems

Integrated power supply and 4 capacitors not only ensure highest, the most stable power output, also being key point to identify true Pico Q-switched YAG Laser (most other competitors have only single capacitor)!!

04 Innovative External Filter

Easy to change and maximize the laser source lifespan!!!

05 Highest energy output

Highest energy output up to 800mj at picosecond mode and 5000mj at long pulse pro mode, ensuring best treatment results!!

06 Interlock

Number one factory with this additional security design, your trustful choice.



LavandulaTM

An Innovative Approach to Best Skin




APPLICATION

- Epidermal Pigmented Lesions
- Enlarge pores
- Acne Scars
- Fine wrinkles
- Uneven Skin Texture / Coloration
- Hair Loss

Technology Behind Lavandula

Lavandula- (Thulium Laser, 1927nm), well known as BB Laser, shows significant brightening outcome.

It is a non-ablative skin resurfacing and rejuvenation laser delivering remarkable results with minimal downtime. Moreover, Lavandula is excellent way to treat Melasma, Pigment lesions, Wrinkles, Large pore and more by adjustable energy intensity and pulse duration simultaneously during procedure. In order to treat pigmentation on Epidermal basal layer, this Sub-ablative from Lavandula system would be the best indication because there are lots of Melanocyte & Melanin on Epidermal basal layer.



Features for Lavandula

- ◆ Highly effective in treatments of skin brightness, whitening, pigmentation, wrinkles, large pore and also improvements in skin texture.
- ◆ Low maintenance cost owing to Non-consumable and long life time of Thulium fiber.
- ◆ Treat diverse indications by separate setting of Pulse duration and Power.
- ◆ Facilitate both characteristics between ablative and non-ablative.

Various sizes of tips and roller tip

Various sizes of tips give users convenience in procedure and roller tip helps better absorption of amplex applying after creation of MTZs (MicroThermal Zones)



MDsculpt iD™

PERSONALIZED SOLUTION

Mdssculpt is the First 6 Bands Monopolar RF and Lipolaser Personalized Body Sculpting Platform that allows you to personalize your ideal body in as little as one 15-minute treatment.



Handpiece

Comprehensive Handpiece Options Maximize Results.



• 6 medium handles (body)



• 3 small handle (face)



• 4 large handle (body)
3 size replacement heads



Android

15.6 Inches Android Control System.



Interface

Intuitive Workstation Interface, 448kHz, 500kHz, 532kHz, 615kHz 680kHz, 1MHz can be freely switched, and cycle work.



Handpieces

13 hand Handpieces, a variety of sizes to meet the needs of different parts.



Temperature

Handle time temperature control, more comfortable and safe.



Experience

Minimal discomfort, no downtime, immediately resume normal activities.



Advantages

Various Applicators with multiple sizes for different treatment areas.

MDSsculpt Technology

MDSsculpt is the first 6 bands monopolar RF and lipolaser personalized body sculpting platform that allows you to personalize your ideal body in as little as one 15-minute treatment. It is a safe and effective technology that is clinically proven to permanently eliminate fat cells in stubborn areas around your abdomen and flanks. With MDSsculpt iD, you can get the body definition you want, even in areas that have been resistant to diet and exercise.



• 3 Small handles (Face)



• 6 Medium handles (Body)



• 4 Large handles (Body)

MDSsculpt Features



Android
15.6 inches android control system.



TCT Control
Innovative Temperature Controlling Technology, more safe and effective.



Intelligent UI
Intuitive Interface with 448kHz, 500kHz, 532kHz, 615kHz, 680kHz, 1MHz working intelligently.



Max Comfort
No pain, no scar and no downtime.



Applicator
Various Applicators with 3 different sizes for different treatment areas.



High ROI
HANDS FREE design and shortest treatment time 15 mins.

RFM Star™

Vacuum Assisted Elevates the
RFM Standard for Skin Solution



APPLICATION


- ◆ Acne scarring
- ◆ Pore reduction
- ◆ Improved skin tone
- ◆ Stretch marks Wrinkle
- ◆ reduction Skin
- ◆ tightening
- ◆ + more



Technology behind RFM STAR

Vacuum Suction Technology for Safe RF Microneedling Treatments

With the application of vacuum suction technology, you can penetrate Radiofrequency energy deep into the dermis to cause a thermal effect. The thermal effect causes an instant contraction to the fibrous septa and dermis which boosts collagen and elastin production, with continuous results.



Unique Vacuum Suction Technology

A clinically proven innovative Vacuum assisted RF micro-needling technology , which increases collagen production in the skin and helps to reduce lines and resurface the skin with virtually no downtime.

With ultra-fine needles and unique vacuum features, it delivers safe and effective Micro-needling treatments with reduced discomfort and downtime compared to other devices.



Before and After



Fractional Micro-Needle RF

It helps you achieve a younger and sleeker look in a lunch time!

Principle:


The machine uses Radio frequency technology to heat the deep, collagen rich layers of your skin. The heat stimulates your body's natural renewal of collagen and promotes the formation of new collagen.

MRF Handle



81pins 49pins 25pins

SRF Handle



SRF

Function

- Face lifting
- Skin tightening
- Wrinkle removal
- Stretch marks removal
- Acne scars removal
- Pore reduction
- Skin texture improvement



Before

After

Before

After

Before

After



EHITTM PLUS

Endogenous heat immunotherapy To get a healthy and sculpted body



- 
Reliable
Temperature Control
- 
Advanced Heating
Method
- 
No Damage to Normal Cells

PRINCIPLE

EHIT Plus is a global leader in the field of Radio-Frequency (RF) solutions for the physiotherapy in aesthetic market. Our ground-breaking technology based on the 448 kHz frequency and Max RF has helped people to get healthy and sculpted body .



R (Resistive Electrode Transfer)

Insulated Electrode
Current to be transferred by Capacity to patient
Go into 3~5cm depth of skin Effective for skin lifting, pain relief



C (Capacitive Electrode Transfer)

Stainless Steel Electrode (not Insulated) Current to be transferred directly to the patient Go into 9~12cm depth of skin Effective for the subcutaneous fat



RF Fat Melting Handle

RF converts electric energy into heat energy to penetrate deep into the skin to accelerate the blood circulation and subcutaneous tissue to stimulate collagen proteins and cell renewal, to tighten skin and slimming effects.





Swan's modelingIndependent Elegant and Small 360°
solid qualitydesktop PortableRemarkable

Principle

Rhea Light mini is light and practical. It is accurate and safe without reducing energy, and has multiple advantages such as "deep penetration, fast effect, short treatment course, lasting effect, and with wide indications" to ensure that the energy of each applicator can evenly penetrate the entire muscle tissue. It has been recognized and widely used by global theaters due to its convenience, non-invasive and quick results.

Configuration advantages

- **Intelligent chip:** precise operation energy is directly without deviation.
- **Advanced temperature sensor:** Intelligent temperature control, to achieve stable output.
- **Golden Electromagnetic Vibration Wave:** Deeply stimulate hormone secretion and enhance quality of life.

Treatment advantages

- **Unique medical applicator**
Strictly selected high-quality medical material, with high energy efficiency on the skin.
- **Smart chip**
The intelligent identification the applicator is new, and high-precision energy output at one time to ensure that the right person for a special use.
- **Safe, hygienic and comfortable**
Medical disposable applicator, safe and hygienic, intelligt temperature control, confort level enhancement.



Effective advantages

- **Immediate effect, Instantly firming**
The clinical effect shows that the firmness is increased by 30% in one treatment, the private health and pelvic floor muscle can be repaired in an all-round way after 3-4 times.
- **Precise and durable recovery**
Remodel the original protein and fascial cells, stimulate the regeneration of nerelet and small blood vessels, it is effectively improving the function of the pelvic fflow with lasting effect.
- **Non-invasive and painless**
There is no downtime and the quality of life enhncement. It is the first choice for thoudands of women for private health.



Immdate effect



Lating improvement



Recovery of
superfical muscle



Recovery of private
health



Recovery and care of
large muscle group